

## ILLUMINATION

### Light Sources

- GS 9200 Halogen light source, 150 watt, focused beam
- GS 9270 Xenon light source, 175 watt, focused beam
- GS 9272 Xenon light source, 300 watt

### Replacement Lamps

- GS 9205 Replacement lamp assembly, 150 watt halogen for Solos ELS-1 (GS 9200)
- GS 9255 Replacement lamp assembly, 175 watt xenon for Solos (GS 9250)
- GS 9256 Replacement lamp assembly, 175 watt xenon for Solos (GS 9480/81) Integrated System
- GS 9259 Replacement lamp assembly, 45 volt metal halide for MP Video light source
- GS 9271 Replacement lamp assembly, 175 watt xenon for Solos GS 9270
- GS 9273 Replacement lamp assembly, 300 watt xenon for Solos GS 9272

### Fiber Optic Cables

- GS 6100 5mm, for interchangeable light source and instrument end (adapters sold separately)
- GS 6100-SS 5mm, Storz light source end to Solos/Storz/Olympus instrument end
- GS 6100-WS 5mm, Wolf light source end to Solos/Storz/Olympus instrument end



*Our Xenon Series is offered in both the 175w and 300w configurations. Fits any fiber optic cable utilizing versatile sleeve adapters. Easy bulb access.*



**Fiber Optic Cables** *(continued)*

---

- GS 6100-OS 5mm, Olympus light source end to Solos/Storz/Olympus instrument end
- GS 6200 5mm, Solos light source end with interchangeable instrument end (adapters sold separately)
- GS 6200-XS 5mm, Solos light source end to Solos/Storz/Olympus instrument end
- GS 6200-XA 5mm, Solos light source end to Circon-ACMI instrument end
- GS 6200-XW 5mm, Solos light source end to Wolf instrument end

NOTE: All in-stock fiber optic cables are 10' length.  
Other lengths available as special order items.

**Light Port Adapters**

---

- GS 9261 For Solos fiber optic cable (black sleeve)
- GS 9262 For Wolf fiber optic cable (gold sleeve)
- GS 9263 For Olympus fiber optic cable (blue sleeve)
- GS 9264 For Storz fiber optic cable (red sleeve)
- GS 9265 For Circon-ACMI/Stryker fiber optic cable (green sleeve)

All light port adapters fit integrated camera/light source combination  
GS 9480, GS 9481, GS 9270, and GS 9271.

**Light Port Adapters (ELS-2 only)**

---


- GS 9251 For Olympus fiber optic cable (blue sleeve)
- GS 9252 For Storz fiber optic cable (red sleeve)
- GS 9253 For Wolf fiber optic cable (gold sleeve)
- GS 9254 For Solos fiber optic cable (black sleeve)
- GS 9257 For Circon-ACMI/Stryker fiber optic cable (green sleeve)

## ILLUMINATION

### Fiber Optic Cable Adapters/Instrument End

---

- GS 4376 For Solos/Storz/Olympus.....
- GS 4377 For Wolf.....
- GS 4378 For Circon-ACMI/Stryker.....



### Fiber Optic Cable Adapters/Light Source End

---

- GS 4380 Male adapter for Storz.....
- GS 4383 Male adapter (long) for Circon-ACMI/Stryker.....
- GS 4386 Male adapter for Wolf/Dyonics.....
- GS 4388 Male adapter for Olympus.....

