

ILLUMINATION

Light Sources

- GS 9200 Halogen light source, 150 watt, focused beam
- GS 9270 Xenon light source, 175 watt, focused beam
- GS 9272 Xenon light source, 300 watt

Replacement Lamps

- GS 9205 Replacement lamp assembly, 150 watt halogen for Solos ELS-1 (GS 9200)
- GS 9255 Replacement lamp assembly, 175 watt xenon for Solos (GS 9250)
- GS 9256 Replacement lamp assembly, 175 watt xenon for Solos (GS 9480/81) Integrated System
- GS 9259 Replacement lamp assembly, 45 volt metal halide for MP Video light source
- GS 9271 Replacement lamp assembly, 175 watt xenon for Solos GS 9270
- GS 9273 Replacement lamp assembly, 300 watt xenon for Solos GS 9272

Fiber Optic Cables

- GS 6100 5mm, for interchangeable light source and instrument end (adapters sold separately)
- GS 6100-SS 5mm, Storz light source end to Solos/Storz/Olympus instrument end
- GS 6100-WS 5mm, Wolf light source end to Solos/Storz/Olympus instrument end



Our Xenon Series is offered in both the 175w and 300w configurations. Fits any fiber optic cable utilizing versatile sleeve adapters. Easy bulb access.



Fiber Optic Cables *(continued)*

GS 6100-OS 5mm, Olympus light source end to Solos/Storz/Olympus instrument end

GS 6200 5mm, Solos light source end with interchangeable instrument end (adapters sold separately)

GS 6200-XS 5mm, Solos light source end to Solos/Storz/Olympus instrument end

GS 6200-XA 5mm, Solos light source end to Circon-ACMI instrument end

GS 6200-XW 5mm, Solos light source end to Wolf instrument end

NOTE: All in-stock fiber optic cables are 10' length.
Other lengths available as special order items.

Light Port Adapters

GS 9261 For Solos fiber optic cable (black sleeve)

GS 9262 For Wolf fiber optic cable (gold sleeve)

GS 9263 For Olympus fiber optic cable (blue sleeve)

GS 9264 For Storz fiber optic cable (red sleeve)

GS 9265 For Circon-ACMI/Stryker fiber optic cable (green sleeve)

All light port adapters fit integrated camera/light source combination
GS 9480, GS 9481, GS 9270, and GS 9271.

Light Port Adapters (ELS-2 only)

GS 9251 For Olympus fiber optic cable (blue sleeve)

GS 9252 For Storz fiber optic cable (red sleeve)

GS 9253 For Wolf fiber optic cable (gold sleeve)


GS 9254 For Solos fiber optic cable (black sleeve)

GS 9257 For Circon-ACMI/Stryker fiber optic cable (green sleeve)

ILLUMINATION

Fiber Optic Cable Adapters/Instrument End

- GS 4376 For Solos/Storz/Olympus.....
- GS 4377 For Wolf.....
- GS 4378 For Circon-ACMI/Stryker.....



Fiber Optic Cable Adapters/Light Source End

- GS 4380 Male adapter for Storz.....
- GS 4383 Male adapter (long) for Circon-ACMI/Stryker.....
- GS 4386 Male adapter for Wolf/Dyonics.....
- GS 4388 Male adapter for Olympus.....

