

# MONITORS

## Sony Medical Monitors

---

PVM 14M2MDU Sony 14" medical monitor,  
600 lines of resolution

PVM 20M2MDU Sony 20" medical monitor,  
600 lines of resolution

*Features:*

- high-resolution Trinitron picture tube
- on-screen menu in five languages
- knob controls of brightness, contrast, aperture, phase, chroma, and volume
- RGB split screen (two RGB inputs)
- analog RGB, Y/C and composite video
- improved contrast CRT
- improved frequency response
- splash-proof panel protects against liquids
- front panel cover prevents inadvertent operation
- automatic selection of standard color systems:  
NTSC and PAL
- UL-2601-1 listed

*Designed specifically for use in medical applications, the M2MD Series monitors provide the exceptional clarity and color fidelity required for endoscopy, ultrasound, surgical and labor microscopy.*