

Video Carts

- GS 9300 Deluxe cart with articulating arm.....
- articulating arm allows up to a 20" monitor
 - adjustable shelves
 - all caster base
front: 5" locking casters; back: 5" casters
 - six-outlet, UL approved, hospital grade electrical strip
 - 12 foot cord on wrap bracket
 - CO₂ tank holder (two tanks)
 - see-through front door
 - locks on front and back doors use same key
 - weight: 295 lbs.
 - overall dimensions: 27"w x 22"d x 60"h



- GS 9370 Endocart with keyboard drawer.....
- accommodates a 14" or 20" monitor
 - adjustable shelves
 - cable management system
 - multi-socket outlet
 - keyboard drawer
 - anti-tamper back panel
 - all caster base
front: 5" locking casters; back: 5" casters
 - weight: 88 lbs.
 - overall dimensions: 25¼"w x 32½"d x 57¼"h



- GS 9380 Mini-cart with keyboard tray.....
- accommodates a 14" or 20" monitor
 - adjustable shelves
 - multi-socket outlet
 - keyboard drawer
 - all caster base
front: 4" locking casters; back: 4" casters
 - 3" locking draw
 - weight: 82 lbs.
 - overall dimensions: 20"w x 24"d x 45"h



MONITOR ACCESSORIES

Video Carts (continued)

- GS 9390 Economy video cart, modular.....
- 5" top module with slots for two velcro straps
 - two 8" modules
 - one locking draw module
 - one 10" module
 - all caster base
 - front: 4" locking casters; back: 4" casters
 - six-outlet, UL approved, hospital grade electrical strip
 - CO₂ tank holder
 - weight: 90 lbs.
 - overall dimensions: 22¼"w x 20"d x 49¼"h

