

## Video Carts

---

- GS 9300 Deluxe cart with articulating arm.....
- articulating arm allows up to a 20" monitor
  - adjustable shelves
  - all caster base
    - front: 5" locking casters; back: 5" casters
  - six-outlet, UL approved, hospital grade electrical strip
  - 12 foot cord on wrap bracket
  - CO<sub>2</sub> tank holder (two tanks)
  - see-through front door
  - locks on front and back doors use same key
  - weight: 295 lbs.
  - overall dimensions: 27"w x 22"d x 60"h



- GS 9370 Endocart with keyboard drawer.....
- accommodates a 14" or 20" monitor
  - adjustable shelves
  - cable management system
  - multi-socket outlet
  - keyboard drawer
  - anti-tamper back panel
  - all caster base
    - front: 5" locking casters; back: 5" casters
  - weight: 88 lbs.
  - overall dimensions: 25¼"w x 32½"d x 57¼"h



- GS 9380 Mini-cart with keyboard tray.....
- accommodates a 14" or 20" monitor
  - adjustable shelves
  - multi-socket outlet
  - keyboard drawer
  - all caster base
    - front: 4" locking casters; back: 4" casters
  - 3" locking draw
  - weight: 82 lbs.
  - overall dimensions: 20"w x 24"d x 45"h



## MONITOR ACCESSORIES

### Video Carts (continued)

---

- GS 9390 Economy video cart, modular.....
- 5" top module with slots for two velcro straps
  - two 8" modules
  - one locking draw module
  - one 10" module
  - all caster base
    - front: 4" locking casters; back: 4" casters
  - six-outlet, UL approved, hospital grade electrical strip
  - CO<sub>2</sub> tank holder
  - weight: 90 lbs.
  - overall dimensions: 22¼"w x 20"d x 49¼"h

